



WHEN YOUR  
ANALYSE SENSITIVITY  
MATTER MOST...




MAYSWAB® Flocked Nasopharyngeal Swab  
Polypropylene, 95 mm Breakpoint

MAYSWAB® Flocked Nasopharyngeal Swab  
Polypropylene, 75 mm Breakpoint

MAYSWAB® Flocked Nasopharyngeal Swab  
Polypropylene, 60 mm Breakpoint

MAYSWAB® Flocked Nasopharyngeal Swab  
Polypropylene, 30 mm Breakpoint



### MAYSWAB® Flocked Nasopharyngeal Swab

It is the ideal product for all rapid diagnostic tests and PCR analysis. Provides excellent absorption capacity with sensitive Flock application. The unique microstructure of flock fibers increases rapid absorption and ensures that biological samples are released at the same speed.

#### QUALITIES

- Perfect Designed Nasopharyngeal Swab
- Superior Absorbent Flock Tip
- Flexible structured Stable Handle
- Optimum use with standard break position
- Made in TURKEY

#### APPLICATIONS

- For Nasopharyngeal sample collection
- Specimen collection
- Rapid diagnostic testing
- Forensic evidence collection



CE ISO 13485



CAT NO	DESCRIPTION
512PP100	MAYSWAB Flocked Nasopharyngeal Swab (Polipropilen, 100 mm Breakpoint)
512PP95	MAYSWAB Flocked Nasopharyngeal Swab (Polipropilen, 95 mm Breakpoint)
512PP80	MAYSWAB Flocked Nasopharyngeal Swab (Polipropilen, 80 mm Breakpoint)
512PP70	MAYSWAB Flocked Nasopharyngeal Swab (Polipropilen, 70 mm Breakpoint)
512PP60	MAYSWAB Flocked Nasopharyngeal Swab (Polipropilen, 60 mm Breakpoint)
512PP50	MAYSWAB Flocked Nasopharyngeal Swab (Polipropilen, 50 mm Breakpoint)
512PP40	MAYSWAB Flocked Nasopharyngeal Swab (Polipropilen, 40 mm Breakpoint)
512PP30	MAYSWAB Flocked Nasopharyngeal Swab (Polipropilen, 30 mm Breakpoint)

MAYSWAB® are produced in accordance with all ISO standards and regulations. We're producing all processes in cleanroom conditions; plastic molding, flocking and also packaging. Products are packaged for protection from contamination and damage in transit for sterile product, individual wrappers meet strict quality standards to assure contents remain sterile as long as the wrapper is intact required regulatory information appears on packaging. This includes LOT number and expiration date, recommended storage temperature, ISO and CE symbols, method of sterilization and additional symbols where appropriate.

ORDERING INFORMATION



Isici Blokleri Mah. Mevlana Bulv. No.182 /B - 84 06530 Cankaya - Ankara - TURKEY  
Phone: 0 312 220 56 50 Fax: 0 312 220 56 51 Email: info@maygen.com.tr